














Item No	Description	Model	Picture	Unit Price
1	<p>EKG</p> <p>ECGs combine all the advanced features of the L-Line models in a compact design. This ECG device is suitable for a broad spectrum of users. Color touch-screen, high-quality ECG printouts and the spirometry option make it a powerful, yet still an easy-to-operate device.</p>	BTC-EC 6		2400
2	<p>Aspirator (Atsos)</p> <p>Portable suction machine for minor surgery, post-surgery treatments, hospital use. Anti-shock case, light, easy to carry and keep clean. Supplied with 1 lt. jar sterilizable, with overflow valve. Complete with suction regulator, vacuum gauge, disposable filter and set of tubes. 230V 50-60Hz.</p>	Fazzini F-20		350
4	<p>Koter</p> <p>400W High Frequency Electrosurgical Unit</p> <p>Technical Specification:</p> <p>Output Frequency: 530 KHzStandard</p> <p>Accessories</p> <ol style="list-style-type: none"> 1.Power cord 2.Single-use ESU pencil with cable 3.Bipolar ESU forceps 4.Single-use Neutral patient plate 5. Double foot switch 6.Reusable bipolar cable 7.Reusable patient plate cable 	DD-2J		3500
5	<p>Table top Pressure Steam Sterilizer</p> <p>Class type: Class B</p> <p>Capacity: 23L</p> <p>Vacuum pump: THOMAS brand</p> <p>Printer: Built-in</p> <p>Spare part: 1 pc door gasket</p> <p>Voltage: 220V, single phase, 50Hz</p>	MST.B(23L)		8560
6	<p>Vertical Pressure Steam Sterilizer</p> <p>Class type: Class B</p> <p>Capacity: 100L</p> <p>Display: Touch screen</p> <p>Printer: Built-in</p> <p>Vacuum pump: Built-in</p> <p>Spare part: 1pc door gasket</p> <p>Voltage: 380V, three phase, 50Hz</p>	Model: MSL.B(100L)		24888

Item No	Description	Model	Picture	Unit Price
7	Vertical Pressure Steam Sterilizer Class type: Class N Model: MSL.N-VS100C Capacity: 100L Display: LED Printer: Without printer Spare part: 1 pc door gasket Warranty: Two years Voltage: 220V, single phase, 50Hz	MSL.N-VS100C		14178
8	Sealing machine Voltage: 220V, single phase, 50Hz			4774
9	Bionet Monitor The Bionet BM5 Patient Monitor is a vital signs monitor with multiple display modes. ECG, SpO2, NIBP, Respiration, 2x IBP* and 2x Temperature* and EtCO2 Ready* 10.4" TFT monitor Individual sweep speed setting for each parameter 128 full hours trend data are stored 20cases * 10 sec alarm waveform are stored, reviewed and printed Li-on rechargeable battery allows cordless use for up to 2 hours	BM5		2300
10	Sigma Amelung Pre packed single macro cuvettes with ball bearing for easy loading Test menu for PT, aPTT, Fibrinogen Store reagent ISI values for automatic INR calculation Single test position to perform all clotting assays Suitable for STAT and routine testing Maintenance free operation	KC 1		800

Item No	Description	Model	Picture	Unit Price
12	<p>Isoflurane</p> <p>The Drager Vapor 2000 Isoflurane is an unheated, calibrated anesthesia vaporizer for the enrichment of dry, fresh medical gas in anesthesia delivery systems with accurately-controlled concentrations of Isoflurane vapor. The Vapor 2000 is connected between the fresh-gas flow-control unit and the fresh-gas outlet.</p> <p>Weight (Filled): 8.8 kg / 19.4 lbs Dimensions (HxWxD): 24.6cm x 10.8cm x 18.8cm Liquid Capacity (Dry Wick): ~360 mL Liquid Capacity (Moist Wick): ~360 mL Flow Range: 0.25 to 15 L/min Maximum Angle of Tilt: 30°</p>	Drager Vapor 2000		1400
14	<p>DocUReader Urine Analyzer</p> <p>This device is both affordable and top-of-the-line. It features a one operation that is simple to use. This device also contains an internal thermal printer. Producing a throughput of up to 40 tests an hour, this device additionally features RS232 serial interface, which makes it useful for hosting either computers or LIS systems.</p>			450
15	<p>Perfusor B. Braun</p> <p>Wide selection of syringes (from 2 - 50 ml) for higher flexibility Infusion rates from 0.01 to 200.0 ml/h Automatic calculation of the infusion rate based on volume over time Stackable with a built-in locking mechanism Battery operable for more than 3 days Automatic recognition of syringe size/fixation Automatic bolus reduction after pressure alarm Occlusion alarm pressure limits selectable: 0.3/0.6/1.2 bar Data Lock Function - keys are blocked</p>	Compact S		800

Item No	Description	Model	Picture	Unit Price
16	<p>Penlon Sigma Delta Vaporizer</p> <p>The Sigma Delta Vaporizer delivers accurate concentrations under varying conditions of flow rate and temperature, particularly at low flow.</p> <p>Anesthetic agent available: Isoflurane or Sevoflurane</p> <p>Filling type available: Key Fill or Funnel Fill (Pour Fill)</p> <p>Mounting type available: Selectatec, Screw Mount or Cagemount</p> <p>Superb performance, particularly at low flow</p> <p>Low body weight</p> <p>Internationally recognized, award –winning design</p>			1700
17	<p>Alma Soprano Ice Laser Diod</p> <ul style="list-style-type: none"> • Virtually painless • Proven Safety Record • All skin types • Tanned skin • Full coverage rate 			28000